



***** YOU MUST VIEW THIS EMAIL FROM YOUR MOBILE DEVICE TO COMPLETE INSTALLATION *****

PowerMic Mobile allows you to use your iPhone or Android smartphone as a secure wireless microphone for dictation into Nuance desktop speech recognition solutions. This application has been developed by Nuance for use with their desktop speech recognition applications. Please follow the instructions below to install, configure and begin using this application.

Step	Instruction	Reference
Step 1:	Download the PowerMic Mobile app from the App Store or Google Play	 iPhone Users (requires iPhone 5 and iOS 15.0 or greater): Android Users (requires Android 4.0 or greater):
Step 2:	Configure the PowerMic Mobile app	Once you have downloaded the app, click the link below from your smartphone to access the PowerMic Mobile configure links (IOS users must copy and paste the link into the browser on the device to launch) IOS: dmic://config/?NmsToken=MjM5MDE0Q0EtMzZCMi00RkQzLUFGNzQtNTUwOEQ2QjgwRTQz Android: http://config/?NmsToken=MjM5MDE0Q0EtMzZCMi00RkQzLUFGNzQtNTUwOEQ2QjgwRTQz
Step 3:	Open PowerMic Mobile	If asked "What application do you want to open the link in," select PowerMic Mobile.
Step 4:	a) In Profile confirm URL b) Profile Name initially appears as Profile 1. You can backspace Profile 1 and apply ROH.	 
Step 5:	PowerMic mobile is used with the following applications.	 DMO Powerchart P710 (For Ambulatory and Inpatient Providers)  DMO FirstNet P710 (For Emergency Medicine Providers)